

BRK™ Transseptal Needles


Product Highlights

- Available in adult and small anatomy sizes
- Bevel angle 50 degrees
- Pointer on the needle shield shows the direction of the needle curve

Ordering Information

Contents: Stainless Steel Transseptal Needle and Stylet (1 unit per box)

Reorder Number	Curve Type	Size	Needle Usable Length (cm)
407200	BRK™	Adult, 18 ga.	71
407201	BRK-1™	Adult, 18 ga.	71
407202	BRK™	Small Anatomy, 19 ga.	56
407203	BRK-2™	Small Anatomy, 19 ga.	56
407205	BRK™	Adult, 18 ga.	89
G407215	BRK-1™	Adult, 18 ga.	89
407206	BRK™	Adult, 18 ga.	98
407207	BRK-1™	Adult, 18 ga.	98



BRK™ Transseptal Needles, XS Series

Product Highlights

- Bevel angle 30 degrees
- Steeper primary bevel angle and two back bevels that combine to form a distinct point at the tip of the needle
- Pointer on the needle shield shows the direction of the needle curve

Ordering Information

Contents: Stainless Steel Transseptal Needle and Stylet (1 unit per box)

Reorder Number	Curve Type	Size	Needle Usable Length (cm)
G407208	BRK™ XS	Adult, 18 ga	71
G407209	BRK-1™ XS	Adult, 18 ga	71
G407210	BRK™ XS	Adult, 18 ga	89
G407216	BRK-1™ XS	Adult, 18 ga	89
G407211	BRK™ XS	Adult, 18 ga	98
G407212	BRK-1™ XS	Adult, 18 ga	98

